



2ml PP vials, 9mm Screw thread ND9



2ml PP vials, 11mm Snap Cap ND11

**Product Features :**

2ml Polypropylene vials are available in 300ul volumes. Breakage is not an issue with PP vials and they are not reactive with most sample chemistries.

- Specification of the bottle is 12 x 32mm, Volume : 300 ul
- Compatible with many solvent.
- Suitable for Agilent, Waters, Brukers, PerkinElmer, ThermoScientifics, Shimadzu, CTC and other rotating or robotic arm samplers.
- Optimal choice when working with proteins.
- Ideal for IC applications.
- Suitable for insert-free microsampling.
- Economical to purchase.
- Translucent and chemically resistant material.
- Traditional the material of choice for LC or LC/MS and CE.
- For pH-sensitive samples, sodium or heavy metals analysis.
- Screw cap and Snap caps available.
- Septa are made of the highest quality materials to ensure proper function



**2ml PP Vials, 9mm Screw Neck ND9**

**Ordering Information :**

**2ml PP Vials, 9mmScrew neck ND9, volume 300ul**

*12 x 32mm, 9mm PP screw thread ND9, Vials*

P/N	Descriptions	Qty
120940	2ml PP screw neck ND9 vial, <b>clear</b> , 300ul	100/pk
120941	2ml PP screw neck ND9 vial, <b>amber</b> , 300ul	100/pk

**9mm Screw Cap and Septa**

P/N	Descriptions	Qty
120907	9mm PP screw thread cap, Blue center hole with Red PTFE/White Silicone septa	100/pk
120908	9mm PP screw thread cap, Blue center hole with White PTFE/Red Silicone septa	100/pk
120909	9mm PP screw thread cap centre hole, 6mm	100/pk
120910	Septa Red PTFE/White Silicone, 9 x 1mm	100/pk
120911	Pre-Slit Septa Red PTFE/White Silicone, 9 x 1mm	100/pk
120912	Septa White PTFE/Red Silicone, 9 x 1mm	100/pk
120914	Pre-Slit Septa White PTFE/Red Silicone pre slit, 9 x 1mm	100/pk
120916	9mm PP screw thread cap, Close-top	100/pk
120917	9mm PP screw thread cap, Blue center hole with Red PTFE/White Silicone/Red PTFE septa	100/pk
120919	9mm PP screw thread cap, Blue center hole with Pre-Slit White PTFE/Red Silicone septa	100/pk
120920	9mm PP screw thread cap, Blue center hole with Pre-Slit Red PTFE/White Silicone septa	100/pk
120924	Septa Red PTFE/White Silicone/Red PTFE, 9 x 1mm	100/pk

**2ml PP vials, 11mm Snap cap ND11**

**Ordering Information:**

**2ml PP Vials, 11mmSnap Cap ND11**

*12 x 32mm, 11mm PP Snap Cap ND11, Vials*

P/N	Descriptions	Qty
121205	2ml PP Snap Cap ND11 vial, <b>clear</b> , 300ul	100/pk
121206	2ml PP Snap Cap ND11 vial, <b>amber</b> , 300ul	100/pk

**PP Snap caps**

*11 mm, Snap Cap 2ml PP Vials, ND11*

P/N	Descriptions	Qty
121228	Natural Seal with White PTFE/Red Silicone Septa	100/pk

**Rack for 2ml Vial**

P/N	Descriptions	Qty
120927	Rack for Vial 2ml, 50 holes	1